Determination of ammonium by spectrophotometer

Baethgen WE, Alley MM (1989) A manual colorimetric procedure for measuring ammonium nitrogen in soil and plant Kjeldahl digests. Commun. Soil Sci. Plant Anal. 20 (9&10) 961-969

Reagents

- 1. Solution of 50% NaOH.
 - a. dissolve 100 g of NaOH in water, dilute to 200 mL
- 2. Buffer solution:
 - a. to a 1000 mL flask add 600 mL of water,
 - b. add 14.2 g of Na2HPO4 or 17.8 g of Na2HPO4.2H2O or 35.8 g of Na2HPO4.12H2O and dissolve
 - c. add 50 g of K Na tartrate (C4H4KNaO6.4H2O)
 - d. add 108 g of 50% NaOH
- 3. Salicylate/nitroprusside solution
 - a. in a 1000mL flask dissolve 150 g of Na Salicylate (C7H5NaO3) and 0.30 g of sodium nitroprusside (Na2[Fe(CN)5NO].2H2O) and make up to 1L
 - b. store in dark in brown bottle
- 4. Hypochlorite solution
 - a. dilute 6 mL of 5.25% sodium hypochlorite to 100 mL
 - b. prepare daily as it isn't stable
- 5. A "Diluent" solution
 - a. This contains the metal catalyst, or digestion solution, or the extraction solution

Standards

- 1. 1000 mg/L NH4-N stock
 - a. In a 1000 mL flask dissolve 4.715 g of ammonium sulphate (NH4)2SO4 (dried at 105C). Make up to 1000 mL with the "diluent" solution
- 2. 100 mg/L NH4-N
 - a. Dilute 10 mL of the 1000 mg/L solution to 100 mL using "diluent"
- 1. Series of NH4 standards
 - a. Dilute 0, 1, 2, 3, 4, 5, 6 of 100 mg/L to 100 mL with "diluent"
 - b. This series contains 0, 1, 2, 3, 4, 5, 6 mg/L of NH4-N

Procedure

- 1. Turn on spectrophotometer and set `wavelength to 650 nm
- 2. Transfer 1mL of solution or standard to a test tube
- 3. Add 5.5 mL of buffer solution. Mix and agitate with vortex
- 4. add 4 mL of salicylate/nitroprusside solution and mix
- 5. add 2 mL of hypochlorite solution and mix
- 6. let rest for 45 min at 25C or 15 min at 37C
- 7. read absorbance at 650 nm within 2 hours